



**TORBAL Class II combination Pharmacy scale and pill counter**  
**- the right solution for every pharmacy**

- NTEP Certified Pill Counting and Compounding.
  - Over 3000 drug database capacity.
- Accurate compounding with readability of 1 mg.
  - Rx and NDC Verification.
  - No cross contamination.
- Meets Handbook 44 requirements for Class A balances.



**TORBAL DRX-4 SERIES**

**Fulcrum Inc.**





## DRX-4C Specifications

Model	DRX-4C
Capacity (Max)	300g
Minimum Load (Min)	0.02g
Reading Unit (d)	0.001g
Verification Unit (e)	0.01g
Tare Range	300g
Accuracy Class	II
Temperature Range	+15C to +30C
Weighing Time	<3s
Pan Dimensions	116mm
Scale Dimensions	235 x 245 x 80mm
Power Supply	Input: 120VDC 60Hz 9.5W Output: 12VDC 500mA
Scale's Net Weight (lb/kg)	8.2 / 3.7
Calibrating Weight	300g



- NDC verification
- Counting transaction receipt printing
- Compounding and recipe-making mode
- Recipe receipt printing
- Handheld Barcode Scanner with automatic activation
- Large LCD Display
- Two RS-232 Interface Ports
- Accurate and easy-to-use leveling system
- Used capacity indicator
- Stabilization and Net Weight indicators
- Center of zero indicator
- Automatic zero-setting mechanism (AZSM)
- Initial zero-setting mechanism (IZSM)
- Sealed front panel
- Draft ring

## Counting Feature Specifications

MW (Minimum piece Weight)	30 mg ( 3 e)
MSS (Minimum Sample Size)	10 pieces
Sample Size Selection	10, 30, 60
Database Capacity	3,000 NDC codes and pill weights



## Legal for trade with NTEP approved counting feature

The DRX-4C is equipped with an NTEP-approved, legal-for-trade counting feature for prescription filling. The scale counts up to 90 pills per second with a database capacity that allows it to store over 3000 NDC codes and pill weights. New drugs and their NDC codes are quickly and easily added to the database, which may be maintained by updating the stored sample sizes and weights.



## NDC verification

The NDC verification feature eliminates potential errors in prescription filling. It ensures that the correct prescription is counted and filled. The feature allows the user to compare two NDC codes and check for a match. The Rx label barcode can be compared to the stock bottle barcode.

## Counting transaction receipt printing

This feature allows the user to print a record of a performed counting transaction. The printed receipt indicates the date and time of the performed transaction as well as the NDC code, average pill weight, and actual count. The printed receipt may be used as a logging tool for the counting feature.

11/23/04 10:49pm

NEW SAMPLE

NAT. DRUG CODE 1191704169

AUG. PILL WT. 0.50[0] g

PILL COUNT 60 PC



TORBAL DRX-4 Series Comparison Chart

Features and Specifications	Model		
	DRX-4C	DRX-4C2	DRX-4
Capacity	300g	300g	100g
Reading Unit (d)	0.001g	0.001g	0.001g
Verification Unit (e)	0.01g	0.01g	0.01g
Tare Range	300g	300g	100g
Accuracy Class	II	II	II
Pan Dimensions	116mm	116mm	116mm
Scale Dimensions	235 x 245 x 80mm	235 x 245 x 80mm	235 x 245 x 80mm
RS-232 Ports	2	1	1
NTEP certified for prescription use	•	•	•
Meets Handbook 44 requirements for a "Class A" prescription Balance	•	•	•
Draft ring	•	•	•
Sealed front panel	•	•	•
Initial zero-setting mechanism (IZSM)	•	•	•
Automatic zero-setting mechanism (AZSM)	•	•	•
Center of zero indicator	•	•	•
Stabilization and NET Weight indicator	•	•	•
Used Capacity indicator	•	•	•
Compounding / Recipe receipt printing	•	•	
NTEP certified and legal for trade Pill Counting feature	•	•	
Counting transaction receipt printing	•	•	
Handheld Barcode Scanner with automatic activation	•		
3,000 drug database capacity	•		
Rx & NDC verification	•		

**ROXON**  
**1-800-361-6991**