

## HBF-306C Fat Loss monitor



- Monitor your progress in losing body fat
- Know you are losing fat, not muscle
- Accurate readings in just seconds
- Athlete mode for accurate results for athletes
- Up to 9 Person Profile

# MODEL HBF-360C

Calculate body fat percentage and Body Mass Index (BMI) in seconds using Bio Impedance Method. Just hold the grip electrodes while standing and press START.

## Measures two Fitness Indicators:

- Body Fat Percentage - Body Mass Index (BMI)

## PRODUCT SPECIFICATIONS...

Model:	HBF-306C OMRON Fat Loss Monitor
Display:	Body fat percentage (4.0 to 50.0%) BI Method BMI (7.0 to 90.0) BMI classification (4 levels)
Set Ranges:	Height: 3' 4" to 6' 6" (100 cm to 199.5 cm) Weight: 23 lbs to 440.5 lbs (10 kg to 199.8 kg) Age: NORMAL; 10 to 80 years old ATHLETE; 18 to 60 years old Gender: Male / Female
Power Supply:	2 AAA batteries (You may use 2 AAA alkaline batteries)
Battery Life:	Approximately 1 year (varies depending on usage and storage.)
Operating Temperature and Humidity:	50°F to 104°F (+10°C to +40°C), 30% to 85% RH
Storage Temperature and Humidity:	-4°F to 140°F (-20°C to +60°C), 10% to 95% RH
External Dimensions:	Approx. 8" (L) x 5" (H) x 2" (W) (197 (L) x 128 (H) x 49 mm (W))
Weight:	Approx. 8 oz. (230 g) (not including batteries)
Components:	1 Fat Loss Monitor, 2 AAA Batteries, 1 Instruction Manual
UPC:	0 73796 30634 2

Specifications and appearance are subject to change without notice for purposes of improvement.

**OMRON**<sup>®</sup>

Distributed by:

**ROXON**  
**1-800-361-6991**

© 3/2006 Omron Healthcare, Inc.

Printed in U.S.A.  
HBF400PS Rev