

# LifeSense®

## WIDESCREE™ Capnography and Pulse Oximetry



### Fast and Easy EtCO<sub>2</sub> and SpO<sub>2</sub> Monitoring

The LifeSense® WIDESCREE™ touch-panel display monitor offers proven PureSAT® SpO<sub>2</sub> and MedAir™ sidestream EtCO<sub>2</sub> technologies for a wide range of clinical settings for both intubated and non-intubated patients. Portable, accurate and cost effective, LifeSense® is perfect for SpO<sub>2</sub> and EtCO<sub>2</sub> monitoring in low acuity areas in and outside the hospital setting.

# LifeSense®

## Product Highlights

- WI DESCREEN™ touch-panel backlit LCD display
- Nonin PureSAT® SpO<sub>2</sub> technology
- MedAir EtCO<sub>2</sub> technology
- Numerical pulse rate and respiratory rate display
- Numerical EtCO<sub>2</sub> and SpO<sub>2</sub> display
- Four-hour trending of pulse rate, SpO<sub>2</sub>, respiration and EtCO<sub>2</sub>
- Capnograph and plethysmograph display
- Audible and visual alarms
- Data output via RS-232 port




## Specifications

Dimensions:	200 (W) x 135 (H) x 50 (D) mm (7.9 x 5.3 x 2 inches)
Weight:	800 grams (1.8 lbs)
<b>Power Data</b>	
Battery charger:	100 – 240 VAC, 50 – 60 Hz
Power consumption:	3.6 W at battery operation, 9 W with power supply
Input:	12 VDC, 720mA
<b>Operation</b>	
Working temperature:	32° to 104°F (0° to +40°C)
Humidity:	10 – 90% (non-condensing)
Atmospheric pressure:	645 – 795 mmHg (860 – 1060 hPa)
<b>Storage</b>	
Storage temperature:	-4° to +122°F (-20° to +50°C)
Humidity:	10 – 90% (non-condensing)
Atmospheric pressure:	Up to 4 atmospheres (83 – 3038 mmHg/110 - 405 hPa)
<b>Pump</b>	
Pump flow:	75 ml/min
Flow accuracy:	±15ml/min
<b>Battery Data</b>	
Type:	Lithium Ion internal battery, rechargeable
Battery capacity:	Approximately 8 hours
Charging time:	Approximately 17 hours or 2 hours for each hour of use
<b>Capnography Measuring Data</b>	
Respiration range:	3 – 80 respirations/min
Update frequency:	Once every breath (No Breath alarm after 25 seconds)
Respiration accuracy:	3 – 50 respirations/min ± 2 51 – 80 respirations/min ± 5
EtCO <sub>2</sub> /CO <sub>2</sub> range:	0 – 99 mmHg, or 0 – 9.9 kPa
EtCO <sub>2</sub> /CO <sub>2</sub> accuracy:	±0.2 kPa / ±2 mmHg, +6% of reading (EtCO <sub>2</sub> /CO <sub>2</sub> reading reaches its steady state accuracy 10 minutes after power up)
Update frequency:	Once every breath (No Breath alarm after 25 seconds)
Trend length:	4 hours; pulse rate, SpO <sub>2</sub> , respiration and EtCO <sub>2</sub>

## Options

- Full line of infant to adult SpO<sub>2</sub> sensors (reusable and disposable)
- Full range of infant to adult CO<sub>2</sub> sampling lines
- Carrying case
- Adjustable mounting clamp
- TrendSense™ memory module, download cable and software
- Digital to analog converter (DAC) cables for pulse and capnography



Infant to Adult CO<sub>2</sub> Sampling Cannulas



T-connectors for Intubated monitoring

## Wide range of PureLight® SpO<sub>2</sub> Sensors

REUSABLE



DISPOSABLE



## Pulse Oximeter Specifications

Displayed Oxygen Saturation Range (SpO <sub>2</sub> )...	0 to 100%
Displayed Pulse Rate Range	18 to 255 beats per minute (BPM)
Measurement Wavelengths***	
Red:	660 nanometers @ 0.8 mW maximum average
Infrared:	910 nanometers @ 1.2 mW maximum average
Saturation Accuracy (Arms) **	70 - 100%

Specifications are subject to change without notice.

Authorized Distributor:

Distributed by:

**ROXON**  
**1-800-361-6991**